











MEDICAL STOOL

9157

INSTRUCTION MANUAL









Technical data	UK	USA
Measurement	Dia 40.5*H54-72cm	Dia 15.9*H21.3-28.3in
Max Load Capacity	150kg	330LB
Angel of seat cushion incline	3°	3°
Forward/backward range of backrest	0~5cm	0~2in
Up and down range of backrest	0~7.5cm	0~3in
Padding	7cm sponge	2.8in sponge




Brush away excess hair



Clean surface with damp cloth or
use simple soap & water solution.



Avoid using corrosive
cleaning agents such as



Do not remove upholstery
for cleaning